**NOT FOR DISTRIBUTION IN THE  
USA OR ITS TERRITORIES**

**NOT295 Rev. D**

## **Diam™ Intervertebral Prosthesis**

ENGLISH

### **Diam™ Intervertebral Prosthesis**

ARTICULAR PROCESS STABILIZATION AND DISCAL ASSISTANCE  
DEVICE

Note: not all parts may be available in each geography.

#### **MATERIAL**

- Dimethyl siloxane (silicone)
- Polyethylene terephthalate
- Titanium

#### **INDICATIONS**

##### **Lumbar Surgery**

- Arthropathic facet-syndrome.
- Posterior prosthesis for discal assistance.
- Transitional osteosynthesis.
- Foraminal stenosis.

The Diam™ Intervertebral Prosthesis allows an effective stabilization of the operated segments while ensuring a certain flexibility. The purpose of this stabilization is to decrease the lumbagos of the patient while avoiding the arthrodesis.

The Diam™ Intervertebral Prosthesis is designed for long-term implantation.

#### **CONTRAINDICATIONS**

Do not implant in case of intolerance to silicone, polyethylene terephthalate, or titanium. Do not implant during child growth.

Do not implant on septic zone or in case of septicaemia.

#### **POSSIBLE UNDESIRABLE SECONDARY EFFECTS**

Known undesirable effects concerning interspinous devices unlikely related to the device, but related to the surgery or patient condition include:

- Inflammation.
- Permanent ligament injury.
- Rupture of ligament.
- Removal of the prosthesis.
- Infection.
- Pseudo meningocele, fistula, persistent CSF leakage, and/or meningitis.
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anaesthesia, paraesthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits.
- Cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss.
- Urinary retention, loss of bladder control, or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles,

and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery.

- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: adjacent segment degeneration.
- Non-union or pseudarthrosis, delayed union.
- Mal-union.
- Cessation of any potential growth of the operated portion of the spine.
- Loss of or increase in spinal mobility or function.
- Inability to perform the activities of daily living.
- Bone loss or decrease in bone density.
- Graft donor site complications including pain, fracture, or wound healing problems.
- Ileus, gastritis, bowel obstruction, or loss of bowel control or other types of gastrointestinal system compromise.
- Haemorrhage, hematoma, occlusion, seroma, oedema, hypertension, embolism, stroke, excessive bleeding.
- Reproduction system compromise, sterility, sexual dysfunction.
- Development of respiratory problems (e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.).
- Change in mental status.
- Death.

Undesirable secondary effects related to the device include:

- Device loosening.
- Device breakage.
- Knitting tearing.
- Seroma/swelling around the device.
- Allergy.
- Foreign body reaction.
- Inflammation.
- Worsening pain.
- Inaccurate positioning (posterior or wrong level).

#### **POSTOPERATIVE**

When explanting and/or disposing of a device, avoid exposure to bodily substances such as blood, tissue, etc., as contact could lead to infection or disease. Always wear and use proper equipment, taking special care with sharp objects and needles. Follow your healthcare center's policy regarding both the disposal of devices and any events of exposure.

#### **STORAGE**

Keep in a dry place, without sunlight, and in an ambient atmosphere.

#### **VISUAL INSPECTION**

Visually inspect all sterile-barrier packaging before use. If the sterile barrier is damaged or the integrity has been compromised, do not use the product. Contact Companion Spine for return information.

Visually inspect the device before use. If the device is damaged, do not use the product. Contact Companion Spine for return information.

#### **PACKAGING**

Devices are supplied sterile. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package is broken, the product should not be re-sterilized. If a loaner set is used, all sets and components should be carefully checked for completeness and to ensure there is no damage prior to use.

Damaged packages or products should not be used, and should be returned to Companion Spine.

#### **WARNINGS**

COUSIN BIOTECH does not guarantee or recommend any special trademark for fixation devices. Properties of these devices are subjected to modifications made by the manufacturer and on which COUSIN BIOTECH has no control.

Diam™ Intervertebral Prosthesis is not recommended for implantation in children. This product must be implanted only by a trained physician.

Physician note: although the physician is the learned intermediary between the company and the patient, the important medical information in this document should be conveyed to the patient. This device was designed for single patient use only. Do not reprocess or reuse this product. Reuse or reprocessing may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

## MRI INFORMATION

In non-clinical testing the **DIAM® Spinal Stabilization System** was determined to be MRConditional. A patient with this device can be safely scanned immediately after device placement under the following conditions:

### **Static Magnetic Field**

- Static magnetic field of 1.5 Tesla and 3-Tesla.
- Maximum spatial gradient magnetic field of 3000-Gauss/cm or less
- Maximum whole body average specific absorption rate (SAR) of 4.0 -W/kg or less under Normal mode or first level controlled operating mode, for 15 minutes of scanning per pulse sequence.

### **MRI-Related Heating**

Based on measurements of RF heating according to ASTM F2182-11a, the DIAM® Spinal Stabilization System is expected to produce a maximum temperature of 3.2°C for a whole body SAR of 4.0 W/kg in a 1.5-Tesla/64-MHz and 3.0-Tesla/128-MHz MRI Systems, for a 15-min of scanning (i.e., per pulse sequence).”

### **Artifact Information**

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the DIAM® Spinal Stabilization System. The maximum artifact beyond the implant was 5 mm for the spin echo sequence and 9 mm for the gradient echo sequence in a 3-Tesla MR system (GE Signa HDxt MR System). Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The presence of other implants or the health state of the patient may require a modification of the MR conditions.

### **MRI Patient Counseling Information**

Physicians should communicate with the patient the following information about Magnetic Resonance Imaging (MRI) with respect to **DIAM® Spinal Stabilization System**:

**DIAM® Spinal Stabilization System** performance has been established for MRI systems at field strengths of 1.5 Tesla and 3.0 Tesla only.

During and MRI the patient may notice a warming sensation around the implant or feel a tingling sensation. If the warming or tingling sensation is uncomfortable the patient should communicate this to the MR technologist, the MRI should be stopped and the settings adjusted to reduce or eliminate the sensation. The highest temperature change above the background observed in non-clinical testing was +3.2°C (associated with specific conditions listed above).

Additionally, the metal in the implant may cause the MRI image to be distorted in the area around the implant. The MRI can be adjusted to minimize the image distortion.

Physicians should instruct patients to:

- Inform any healthcare personnel (e.g. doctor or MR technologist) that they have an implanted DIAM® Spinal Stabilization System prior to receiving an MRI.
- The patient’s doctor will recommend whether or not an MRI is appropriate.

### **PRODUCT COMPLAINTS**

To report product problems, contact Companion Spine.

Patients in the European Union experiencing a serious incident in relation to the device should contact Companion Spine and the competent authority of the Member State in which they are established.

## FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is required, contact Companion Spine.  
©2022 Companion Spine. All rights reserved.

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