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**NOT068 / 240531**

Version du 31/05/2024

## **REINFORCEMENT ARTICULAR LIGAMENTS STERILE ANCHOR DEVICES STERILE PRODUCT FOR SINGLE USE**

### **DESCRIPTION**

Reinforcing Ligaments for orthopedic surgery.

### **MATERIALS**

Polyethylene terephthalate (No resorbable)

Stainless steel

Nor human neither animal origin.

### **INDICATIONS**

The ligaments are intended for the surgical reinforcement of instabilities of the various joints. They can be implanted in intra or extra-articular. The reinforcement ligament must be always used in addition to suture or natural ligament graft.

### **CONTRAINDICATIONS**

Do not use:

- In a septic environment or in an infected site,
- allergy to one of components
- growing child.

### **POSSIBLE ADVERSE SIDE EFFECTS**

- Routine surgical complications
- Superficial or deep infection
- Anesthesia complications
- Local hematoma, blood loss
- Postoperative loss of motion, frozen joint
- Inflammatory / allergic reactions / arthro fibrosis (around 1%) / Bone resorption (osteolysis)
- Ongoing instability episodes
- Neurovascular injury, peripheral neuropraxia
- Wound dehiscence, serous leak
- Post-surgical pain, persistent pain
- Ligament injury
- Permanent elongation of the ligament / Rerupture of the ligament

### **PRECAUTIONS FOR USE**

Please check the perfect integrity of packaging (do not use if the labels and/or the protective packaging are damaged). Do not use if the device is damaged.

**DO NOT REUSE**

**DO NOT RESTERILIZE**

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these devices are the subject of modifications by the manufacturer and over which COUSIN BIOTECH has no control. Suture clips or sharp mounting screws must not be used.

This device should be implanted **only** by a qualified surgeon who has been trained in the use of this product.

#### *Postoperative precautions*

The device has not been tested for heating or migration in the MR environment.

The device has been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment.

According to medical literature and in view of the materials used in the implant, the device can be considered as MR safe.

### **PRECAUTIONS FOR STORAGE**

To be stored in a dry place, protected from sunlight, at room temperature in the original packaging.

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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