ADJUSTABLE ACL CORTICAL FIXATION

STERILE CORTICAL FIXATION DEVICES

en	<u>Instructions for use</u>	page	2
fr	<u>Notices d'instructions</u>	page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
pt	<u>Nota de instruções</u>	Pagina	12
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	14
bg	<u>Инструкции за употреба</u>	Страница	16
tr	<u>Kullanım Kilavuzu</u>	Sayfa	19
fi	<u>Käyttöohjeet</u>	Sivu	21
ru	<u>Инструкция по применению</u>	Страница	23

Distributed by:

bioretec

BIORETEC Ltd.
Yrittäjänkulma 5,
FI-33710 Tampere, Finland
Tel. : +358 20 778 9500
www.bioretec.com

COUSIN BIOTECH
Allée des Roses
59117 Wervicq-Sud – France
Tel. : +33 (0) 3 20 14 41 20
Fax : +33 (0) 3 20 14 40 13
www.cousin-biotech.com

CE
1639

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**ADJUSTABLE ACL CORTICAL FIXATION
STERILE CORTICAL FIXATION DEVICES
SINGLE USE STERILE PRODUCT**

DESCRIPTION

Cortical Fixation System for tendons and ligaments.

The adjustable fixation System comprises an adjustable loop for graft suspension and a titanium button for the cortical support.

MATERIALS

- *Materials implanted:*

Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Adjustable loop: 100% polyethylene (UHMWPE)

Extended plate for cortical fixation System (optional): Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Neither animal nor human origin. Non-absorbable

- *Materials in contact with the patient during the operation:*

Traction threads: polyethylene terephthalate + biocompatible colorant D&C green#6

INDICATIONS

The Cortical Fixation Systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures for the anterior cruciate ligament (ACL) reconstruction.

PERFORMANCE

This SYSTEM facilitates ligament or tendon repair/reconstruction by graft fixation in the bone.

It has been approved for long-term implantation in the human body.

CONTRAINDICATIONS

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply and previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Not to be used for the bone – patella tendon – bone technique.
- Do not implant in children while they are growing.

POSSIBLE ADVERSE EFFECTS

Possible complications linked to the implant:

- 1- Breaking of the loop, breaking of the traction wires, breaking of the button
- 2- Failure when placed under tension
- 3- Migration or upturn of the button
- 4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

- 1- Pain in the incision area or surgical site
- 2- Local and transient inflammatory reactions
- 3- Infection
- 4- Wound dehiscence

WARNING

- 1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.
- 2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the SYSTEM before the bone has sufficiently healed.
- 3- Ensure that the bone tunnel has the correct length and diameter.

PRECAUTIONS BEFORE USE

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the SYSTEM and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN BIOTECH recommends that instruments should be regularly inspected for wear and bending.

DO NOT REUSE

DO NOT RESTERILIZE

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN BIOTECH does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN BIOTECH has no control.

This device must **only** be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

SURGICAL TECHNIQUE

Preliminary information

The system is adapted for the use of the right internal and semitendinosus tendons with the all-inside and inside-out techniques as well as retrodrilling.

In case of retrodrilling, COUSIN BIOTECH recommends a tunnel diameter of 4.5 mm to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 4.5 mm has not been validated. If the tunnel diameter is greater than 4.5 mm, the surgeon can decide whether or not to use the extended plate.

Technique

1. Suture the graft and put it in the loop.
2. Create a bone tunnel to insert the implant with a diameter of 4.5 mm (recommended diameter) and a graft cavity suited to its diameter.
3. Insert the suture threads and adjustment threads through the bone tunnel.
4. Pull on one of the sutures until the button comes out of the bone tunnel through the cortical bone whilst maintaining the suture and adjustment threads.
5. Pull on the second suture in order to tip over the button and allow cortical contact.
6. Validate fixation by pulling the graft in the opposite direction.
7. While still applying counter tension to the graft, slide the adjustment threads one after another to adjust the loop to the required length, starting with the white thread. Do not apply tension to the black thread.
8. To lock the loop, continue to apply tension to the white adjustment thread and pull firmly on the black adjustment thread.
9. Fix the graft in the opposite bone tunnel using the required method.
10. If necessary, slide the adjustment threads alternately once again to obtain the required tension. And repeat the step 8.
11. Remove the suture thread(s) from the button.
12. Section the adjustment threads at a minimum distance of 10 mm from the cortical button.

Use of the extended plate

It may be necessary to use an extended plate to increase the contact surface of the button on the cortical bone (full tibial or femoral tunnel). Repeat steps 1 to 5 above.

6. Remove one of the two sutures.
7. Put the button on the extended plate. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.
8. While still applying counter tension to the graft, slide the adjustment threads one after another to adjust the loop to the required length, starting with the white thread. Do not apply tension to the black thread.
9. To lock the loop, continue to apply tension to the white adjustment thread and pull firmly on the black adjustment thread.
10. Fix the graft in the opposite bone tunnel using the required method.
11. If necessary, slide the adjustment threads alternately once again to obtain the required tension. Repeat step 9.
12. Remove the remaining suture thread from the button.
13. Section the adjustment threads at a minimum distance of 10 mm from the cortical button.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight, at room temperature and in the original packaging.

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

IMPLANT CARD

Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.

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