COMETE ®

FIXATION FEMORALE CORTICALE LCA

Instructions for use

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Version du 03/04/2024

Rx
only

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician

This release is the last update of the instructions of use and replace the previous edition

COMETE[®]
ACL FEMORAL CORTICAL ATTACHMENT
STERILE CORTICAL FIXATION DEVICES
SINGLE USE STERILE PRODUCT

DESCRIPTION

Cortical Fixation System for tendons and ligaments.
The cortical fixation system comprises a continuous loop for graft suspension and a titanium button for the cortical support.

MATERIALS

- *Materials implanted:*

Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI)
Continuous loop: 100% polyethylene terephthalate (PET)
Extended plate for cortical fixation system (optional): Non-ferromagnetic titanium alloy (Ti TA6V ELI)
Neither animal nor human origin. Non-absorbable

- *Materials in contact with the patient during the operation:*

Traction threads: polyethylene terephthalate + biocompatible colorant D&C green#6

INDICATIONS

The Cortical Fixation Systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures for the anterior cruciate ligament (ACL) reconstruction.

PERFORMANCE

This system facilitates ligament or tendon repair/reconstruction by graft fixation in the bone.
It has been approved for long-term implantation in the human body.

CONTRAINDICATIONS

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply and previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Not to be used for the bone – patella tendon – bone technique.
- Do not implant in children while they are growing.

POSSIBLE ADVERSE EFFECTS

Possible complications linked to the implant:

- 1- Breaking of the loop, breaking of the traction wires, breaking of the button
- 2- Failure when placed under tension
- 3- Migration or upturn of the button
- 4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

- 1- Pain in the incision area or surgical site
- 2- Local and transient inflammatory reactions
- 3- Infection
- 4- Wound dehiscence

WARNING

- 1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.
- 2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the system before the bone has sufficiently healed.
- 3- Ensure that the bone tunnel has the correct length and diameter.

PRECAUTIONS BEFORE USE

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the system and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN BIOTECH recommends that instruments should be regularly inspected for wear and bending.

DO NOT REUSE

DO NOT RESTERILIZE

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN BIOTECH does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN BIOTECH has no control. This device must **only** be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

SURGICAL TECHNIQUE

Preliminary information

The system is adapted for the use of the right internal and semitendinosus tendons with the all-inside and inside-out techniques as well as retrodrilling.

In case of retrodrilling, COUSIN BIOTECH recommends a tunnel diameter of 5 mm to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 5 mm has not been validated. If the tunnel diameter is greater than 5 mm, the surgeon can decide whether or not to use the extended plate.

Technique

1. Suture the graft and put it in the loop.
2. Create a bone tunnel to insert the implant with a diameter of 5 mm (recommended diameter) and a graft cavity suited to its diameter.
3. Choosing the size of the implant. Total length of the femoral tunnel (between the two cortices) - length of the graft tunnel +5 cm. If the result is not equal to the size of the implant, round up to the larger implant size.
4. Pass the threads through the bone tunnel.
5. Pull on one of the sutures until the button comes out of the bone tunnel through the cortical bone.
6. Tip over the button against the cortical bone surface by pulling on the graft.
7. Fix the graft in the opposite bone tunnel using the required method.
8. Remove the suture thread(s) from the button.

Use of the extended plate

It may be necessary to use an extended plate to increase the contact surface of the button on the cortical bone (full tibial or femoral tunnel). Repeat steps 1 to 6 above.

7. Remove one of the two sutures.
8. Slide the extended plate under the button. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.
9. Fix the graft in the opposite bone tunnel using the required method.
10. Remove the remaining suture thread from the button.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight, at room temperature and in the original packaging.

EXPLANTATION AND DISPOSAL OF THE DEVICE

COUSIN BIOTECH recommends that explantation and handling should be carried out in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

All explanted devices must be returned for analysis according to the protocol in force. This is available on request from COUSIN BIOTECH. It should be noted that any implants that cannot be cleaned and disinfected before shipping must be placed in a sealed package.

Disposal of explanted medical devices must be carried out in accordance with the standards in force in the country, concerning the disposal of infectious waste.

There are no special recommendations for the disposal of non-implanted devices.

INFORMATION REQUEST AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH undertakes to make every effort to produce and supply a high quality medical device. However, if a health professional (customer, user, consultant etc.) has a complaint or grounds for dissatisfaction with a product, regarding its quality, safety or performance, he/she must inform COUSIN BIOTECH as soon as possible.

In the event of malfunctioning of an implant or if it contributes towards a serious adverse effect for the patient, the healthcare centre must follow the legal procedures in force in the country, and immediately inform COUSIN BIOTECH.

In all correspondence, please state the reference, batch number, name and address of a consultant, and provide a full description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

COMETE[®]
FIXATION FEMORALE CORTICALE LCA
DISPOSITIFS DE FIXATION CORTICALE STÉRILES
PRODUIT STERILE A USAGE UNIQUE

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Patient information leaflet COMETE[®]

Item	Patient information leaflet
1	a. The name of the device : COMETE [®] b. The model of the device : OAMGEFX15U, OAMGEFX20U, OAMGEFX25U, OAMGEFX30XU, OAMGEFX35U, OAMGEFX40U
2	a. The intended purpose of the device : The Cortical Fixation Systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures for the anterior cruciate ligament (ACL) reconstruction. b. The kind of patient on whom the device is intended to be used : Patient with Anterior Cruciate Ligament (ACL) reconstruction
3	Special operating instructions for the use of the device : The device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
4	a. The intended performance of the device : This system facilitates ligament or tendon repair/reconstruction by graft fixation in the bone. It has been approved for long-term implantation in the human body. b. Undesirable side effects : Possible complications linked to the implant: <ul style="list-style-type: none"> • Breaking of the loop, breaking of the traction wires, breaking of the button • Failure when placed under tension • Migration or upturn of the button • Reaction to a foreign body Possible complications linked to insertion of the implant: <ul style="list-style-type: none"> • Pain in the incision area or surgical site • Local and transient inflammatory reactions • Infection • Wound dehiscence
5	When COMETE[®] Should Not be Used : <ul style="list-style-type: none"> • Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.

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5	<ul style="list-style-type: none"> • Insufficient quantity or quality of bone. • Blood supply and previous infections which may tend to retard healing. • Active infection. • Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period. <ul style="list-style-type: none"> • Not to be used for the bone – patella tendon – bone technique. • Do not implant in children while they are growing.
6	<p>MRI Information : The implanted raw materials used are non-ferromagnetic materials.</p>
7	<p>Precautions :</p> <p>1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.</p> <p>2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the system before the bone has sufficiently healed.</p> <p>3- Ensure that the bone tunnel has the correct length and diameter.</p>
8	<p>Materials included in the device:</p> <ul style="list-style-type: none"> - Materials implanted: Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI) Continuous loop: 100% polyethylene terephthalate (PET) Extended plate for cortical fixation system (optional): Non-ferromagnetic titanium alloy (Ti TA6V ELI) Neither animal nor human origin. Non-absorbable - Materials in contact with the patient during the operation: Traction threads: polyethylene terephthalate + biocompatible colorant D&C green#6
9	<p>Reporting a serious incident :</p> <p>If you encounter any serious incident with your COMETE[®] implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov</p>
10	<p>Legal Manufacturer and Manufacturing Facility : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille