# PROMESH® SURG ABSO VENT

## STERILE SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	4
de	<a href="#">Gebrauchsanweisung</a>	Seite	6
it	<a href="#">Istruzioni per l'uso</a>	Pagina	8
es	<a href="#">Instrucciones de uso</a>	Pagina	10
pt	<a href="#">Instruções de utilização</a>	Pagina	12
ελ	<a href="#">Οδηγίες χρήσης</a>	Side	14
cs	<a href="#">Návod k použití</a>	Strana	16
da	<a href="#">Anvendelsesinstruktioner</a>	Side	18
fi	<a href="#">Käyttöohje</a>	Sivu	20
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	22
pl	<a href="#">Instrukcja użytkowania</a>	Strona	24
ro	<a href="#">Instrucțiuni de utilizare</a>	Pagină	26
sk	<a href="#">Návod na použitie</a>	Strana	28
sl	<a href="#">Navodila za uporabo</a>	Stran	30
tu	<a href="#">Kullanım talimatları</a>	Sayfa	32
sv	<a href="#">Bruksanvisning</a>	Sida	34
hu	<a href="#">Használati utasítás</a>	Oldal	36
no	<a href="#">Bruksanvisning</a>	Side	38
ru	<a href="#">Инструкция по применению</a>	страница	40
uk	<a href="#">Інструкція із застосування</a>	сторінка	42
cn	<a href="#">使用说明</a>	页	44
ar	<a href="#">تعليمات الاستخدام</a>	صفحة	46

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This release is the last update of the instructions of use and replace the previous edition

**PROMESH® SURG ABSO VENT**

STERILE SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
**EXTRAPERITONEAL**  
SINGLE USE STERILE PRODUCT

**DESCRIPTION**

PROMESH® SURG ABSO VENT meshes are extraperitoneal partially resorbable parietal reinforcement implants.

**IMPLANTED MATERIALS**

Polypropylene (non-resorbable knitted monofilament) – Poly L Lactic Acid (resorbable knit)  
Origin neither human nor animal – Semi resorbable.

**INDICATIONS**

Repair of ventral hernia or other fascial defects that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

**PERFORMANCE**

The PROMESH® SURG ABSO VENT implant is designed for extraperitoneal implantation only. Meshes are biocompatible and partially resorbable. They have the advantage of being cut-able, of having a very high resistance to suture, of having a great flexibility.

**CONTRAINDICATIONS**

Do not implant intraperitoneally.

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Recurrence
- Hematoma
- Foreign body reaction
- Infection
- Adhesion formation
- Mesh migration
- Irritation nearby organ
- Inflammation
- Erosion
- Mesh deformation
- Lymphocele
- Fistula formation
- Seroma
- Allergic reaction

**USAGE PRECAUTIONS**

PROMESH® SURG ABSO VENT implants are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labeling, the PROMESH® SURG ABSO VENT implants are for single use only. They can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

**EXPLANTATION AND ELIMINATION OF DEVICES**

Explanations and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected before dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local

regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUEST AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPORTANT**

For more information about the use of this product, please contact your representative or distributor.


### **SURGICAL PROCEDURE**

#### **Extra-Peritoneal Hernia repair**

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents.

- Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.
- Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath.
- The **PROMESH® SURG ABSO VENT** implant is placed in the newly formed pre-peritoneal space. **PROMESH® SURG ABSO VENT** mesh has an oriented elasticity:

Direction of the elasticity

- Fixation of the mesh if needed, to the muscle layer above
- The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents.
- Redon's drains under aspiration are usually placed. The anterior sheath is closed. If tension is present, relaxing incisions may be used
- Skin closure.
- Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour.
- A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

**R<sub>x</sub> only** Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

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