PROMESH® SURG UMB LIGHT

STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
pt	Instruções de utilização	Pagina	12
nl	Gebruiksaanwijzing	Pagina	14
ελ	Οδηγίες χρήσης	Side	16
ro	Instrucțiuni de utilizare	Pagină	18
pl	Instrukcja użytkowania	Strona	20
sl	Navodila za uporabo	Stran	22
cs	Návod k použití	Strana	24
da	Anvendelsesinstruktioner	Side	26
fi	Käyttöohje	Sivu	28
tr	Kullanım talimatları	Sayfa	30
sv	Bruksanvisning	Sida	32
sk	Návod na použitie	Strana	34
hu	Használati utasítás	Oldal	36
no	Bruksanvisning	Side	38
ru	Инструкция по применению	страница	40
uk	Інструкція із застосування	сторінка	42
cn	使用说明	页	44
ar	تعليمات الاستخدام	صفحة	46

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This release is the last update of the instructions of use and replace the previous edition

PROMESH® SURG UMB LIGHT

STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

Expandable parietal prosthesis with inflatable balloon for visceral surgery

STERILE PRODUCT FOR SINGLE USE

EXTRA-PERITONEAL USE

DESCRIPTION

The **PROMESH® SURG UMB LIGHT** prosthesis consists of:

- A parietal prosthesis made of a plug in polypropylene and two polyethylene terephthalate and ePTFE fixation threads crimped with stainless steel needles.
- A balloon which allows the deployment of the parietal prosthesis (by using the syringe included in the medical device).

IMPLANTED MATERIALS

Polypropylene, Polyethylene terephthalate, Biocompatible dye (DC Green #6), ePTFE,

MATERIALS IN CONTACT WITH THE PATIENT DURING SURGERY:

Stainless steel, Poly vinyl chloride.

Origin neither human nor animal – Non-absorbable.

INDICATIONS

The **PROMESH® SURG UMB LIGHT** is intended for the use in all forms of hernia repair requiring reinforcement with non-absorbable support material

PERFORMANCES

The **PROMESH® SURG UMB LIGHT** medical device is particularly designed to be extra-peritoneally implanted for the mini invasive surgical treatment (small incision) of all form of hernia.

The mesh is biocompatible and non-absorbable. It has the advantage of having a very high resistance to suturing, and a great flexibility.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant could generate possible undesirable effects:

- | | | |
|---------------------|---------------------------|------------------|
| - Discomfort/Pain | - Recurrence | - Adhesions |
| - Erosions | - irritation nearby organ | - Infection |
| - Inflammation | - Fistula formation | - Seroma |
| - Mesh deformation | - Hematoma | - Mesh migration |
| - Allergic reaction | - Foreign body reaction | - Lymphocele |

PRECAUTIONS FOR USE

PROMESH® SURG UMB LIGHT devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

IT IS IMPERATIVE THAT THE BALLOON IS WITHDRAWN


As indicated on the product labelling, the **PROMESH® SURG UMB LIGHT** implant is designed for single use. It should in no event be reused and/or reesterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

STORAGE OF THE PROSTHESIS

Store in a dry place away from sunlight and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

SURGICAL PROCEDURE

SURGICAL PROCEDURE				
<p>1-Installation: - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro muscular</p>		<p>2-Deployment: - The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe. - Parietal fixing is carried out by two points thanks to the pre-positionned sutures.</p>		<p>3- Prosthesis in place: -The prosthesis is fixed, -This balloon must imperatively be withdrawn</p> 

SURGICAL TECHNIQUE FOR UMBILICAL HERNIA REPAIR

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space.. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis.

Be careful not to suture the balloon to the implant when fixing the prosthesis.

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

EXPLANTATION AND ELIMINATION OF DEVICES

Explanations and handling should be carried out in accordance with ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected before dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUEST AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

Caution : Federal law (USA) restricts this device to sale,distribution and use by or on the order of a physician.

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