Ligament Acromio-claviculaire Acromioclavicular ligament Ligamento Acromioclavicular

ACROLIG

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01440 Viriat – France

Tél. : (33) 04 74 55 35 55 – Fax : (33) 04 74 52 44 01

E-mail : info@fxsolutions.fr – www.fxsolutions.fr

COUSIN BIOTECH

Allée des Roses
59117 Wervicq-Sud -
FRANCE

Tél. : +33 (0) 3 20 14 41 20

Fax : +33 (0) 3 20 14 40 13

www.cousin-biotech.com

Made in France

CE
1639

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NOT198 / 240711

Version du 11/07/2024

ACROLIG ACROMIOCLAVICULAR LIGAMENT NON RESORBABLE ARTICULAR LIGAMENTS STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

Reinforcing Ligaments for orthopedic surgery.

MATERIALS

Polyethylene terephthalate– Titanium TA6V ELI - ISO 5832-3

Nor human neither animal origin.

Non resorbable.

INDICATIONS

The ligaments are intended for the surgical reinforcement of instabilities in different joints. They can be located intra- or extra-articularity, by arthrotomy or arthroscopy. The textile structure enables a good recolonization of the tissues. The reinforcing ligaments should always be used to complement a suture or natural ligament graft.

CONTRAINDICATIONS

Do not use:

- In a septic environment or in an infected site,
- allergy to one of components
- growing child.

POSSIBLE ADVERSE SIDE EFFECTS

- Inflammatory reactions
- Permanent pulled ligament
- Rupture of the product
- Clavicular bony erosion
- Failure of the coracoid /clavicular button
- Fracture of the distal clavicle at the clavicular hole
- Arthritis of the acromioclavicular joint

PRECAUTIONS FOR USE

Please check the perfect integrity of packaging (do not use if the labels and/or the protective packaging are damaged).

Do not use if the device is damaged.

DO NOT REUSE

DO NOT RESTERILIZE

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these devices are the subject of modifications by the manufacturer and over which COUSIN BIOTECH has no control. Suture clips or sharp mounting screws must not be used.

This device should be implanted **only** by a qualified surgeon who has been trained in the use of this product.

Postoperative precautions

The device ACROLIG has not been tested for heating or migration in the MR environment.

The device has been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment.

According to medical literature and in view of the materials used in the implant, the device can be considered as MR safe.

PRECAUTIONS FOR STORAGE

To be stored in a dry place, protected from sunlight, at room temperature in the original packaging.

EXPLANATION AND ELIMINATION OF DEVICES

Explanations and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected before dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUEST AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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