en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notices d'instructions</a>	Page	4
de	<a href="#">Gebrauchsanweisung</a>	Seite	6
it	<a href="#">Istruzioni per l'uso</a>	Pagina	8
es	<a href="#">Instrucciones de uso</a>	Pagina	10
pt	<a href="#">Nota de instruções</a>	Pagina	12
ελ	<a href="#">Οδηγίες χρήσης</a>	σελίδα	14
ro	<a href="#">Instrucțiuni de utilizare</a>	Pagină	16
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	18
da	<a href="#">Brugsvejledning</a>	Side	20
fi	<a href="#">Käyttöohjeet</a>	Sivu	22
hu	<a href="#">Felhasználási útmutató</a>	Pagina	24
pl	<a href="#">Instrukcja użycia</a>	Strona	26
sk	<a href="#">Návod k použiti</a>	Strana	28
cs	<a href="#">Pokyny k použiti</a>	Strana	30
tu	<a href="#">Kullanım Kilavuzu</a>	Sayfa	32
ru	<a href="#">Инструкция по применению</a>	страница	34
cn	<a href="#">使用说明</a>	页	36
ar	<a href="#">تعليمات الاستخدام</a>	صفحة	38
sl	<a href="#">Navodila za uporabo</a>	Stran	39
sv	<a href="#">Bruksanvisning</a>	Sida	41
no	<a href="#">Bruksanvisning</a>	Side	43
uk	<a href="#">Інструкція із застосування</a>	сторінка	45

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**NOT281\_240917**  
 Version du 17/09/2024

This release is the last update of the instructions of use and replace the previous edition

**PROMESH® SURG FLEXIL**  
**Sterile extraperitoneal non-resorbable parietal reinforcement implant**  
**STERILE PRODUCT FOR SINGLE USE**

**DESCRIPTION**

**PROMESH® SURG FLEXIL** meshes are parietal reinforcement implants.

**MATERIALS**

Polyethylene terephthalate (see on the package)

Polyethylene terephthalate deep-dyed with dimethyl siloxane (see on the package)

Marked thread: Polyethylene terephthalate tinted with a biocompatible dye (depending on references) (See on the package).

**Origin neither human nor animal – Non-absorbable.**

**INDICATIONS**

Repair and parietal reinforcement for inguinal, femoral and ventral hernias.

**PERFORMANCES**

Meshes are adapted to both laparoscopic and conventional surgery.

The **PROMESH® SURG FLEXIL** device is designed for extraperitoneal implantation.

The mesh is biocompatible and non-absorbable, causing a fibrosis reaction which takes over the reinforcement after six months of implantation. It has the advantage of having a shape memory, being cut-able, having a very high resistance to suturing, and great flexibility.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy
- Surgery via a vaginal approach

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- Discomfort/Pain
- Adhesion formation
- Hematoma
- Foreign Body Reaction
- Infection
- Fistula formation
- Mesh migration
- Irritation near by organ
- Recurrence
- Erosion
- Mesh deformation
- Inflammation
- Seroma
- Allergic reaction

**PRECAUTIONS FOR USE**

**PROMESH® SURG FLEXIL** meshes are delivered sterile.

Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN SURGERY does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a **qualified** surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE**

As indicated on the product labelling, the implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

### **STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

**For more information about the use of this product, please contact your representative or distributor.**

### **[Table of contents](#)**

**R<sub>X</sub>**  
only

**Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.**

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