en	Instructions for use	Page	2
fr	Notice d'instructions	Page	5
de	Gebrauchsanweisung	Seite	8
it	Istruzioni per l'uso	Pagina	11
es	Instrucciones de uso	Pagina	14
pt	Instruções de utilização	Pagina	17
ελ	Οδηγίες χρήσης	Side	20
da	Anvendelsesinstruktioner	Side	23
fi	Käyttöohje	Sivu	25
nl	Gebruiksaanwijzing	Pagina	27
ro	Instrucțiuni de utilizare	Pagină	30
tr	Kullanım talimatları	Sayfa	33
pl	Instrukcja użytkowania	Strona	35
sl	Navodila za uporabo	Stran	38
cs	Návod k použití	Strana	40
sv	Bruksanvisning	Sida	43
sk	Návod na použitie	Strana	45
hu	Használati utasítás	Oldal	48
no	Bruksanvisning	Side	50
ru	Инструкция по применению	страница	52
uk	Інструкція із застосування	сторінка	55
zn	使用说明	页	58
ar	تعليمات الاستخدام	صفحة	61

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NOT284_240918
Version du 18/09/2024

This release is the last update of the instructions of use and replace the previous edition

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

PROMESH® SURG UMB

EXPANDABLE PARIETAL PROSTHESIS WITH INFLATABLE BALLOON FOR VISCERAL SURGERY SINGLE USE STERILE PRODUCT

DESCRIPTION

The **PROMESH® SURG UMB** prosthesis consists of:

- A parietal prosthesis made of an ePTFE and polypropylene plug, two or four polyethylene terephthalate with biocompatible colorant (D&C Green #6) and ePTFE fixation threads crimped with stainless steel needles (see on package).
- A balloon which allows the deployment of the parietal prosthesis (except for PSUMBIN05D).

IMPLANTED MATERIALS :

Polypropylene, Polyethylene terephthalate, ePTFE, Biocompatible colorant (D&C Green #6).

MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY :

Stainless steel, Poly vinyl chloride.

Origin neither human nor animal – Non-absorbable.

INDICATIONS:

Repair and parietal reinforcement for umbilical and ventral hernias.

PERFORMANCES

The **PROMESH® SURG UMB** prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for the mini invasive surgical treatment of umbilical and ventral hernias.

Products, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow a optimum and rapid integration and colonisation.

If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- | | | |
|---------------------|---------------------------|------------------|
| - Discomfort/Pain | - Recurrence | - Adhesions |
| - Erosions | - irritation nearby organ | - Infection |
| - Inflammation | - Fistula formation | - Seroma |
| - Mesh deformation | - Hematoma | - Mesh migration |
| - Allergic reaction | - Foreign body reaction | - Lymphocele |

PRECAUTIONS FOR USE

PROMESH® SURG UMB prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN SURGERY does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE THIS BALLOON MUST IMPERATIVELY BE WITHDRAWN

As specified on the product's labeling, the **PROMESH® SURG UMB** mesh is for single use only. It cannot be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your representative or distributor.

SURGICAL PROCEDURE (except for PSUMBIN05D which comes without balloon)					
1-Installation: - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro parietal (intraperitoneal)		2-Deployment: - The expansion and the spreading of the prosthesis are obtained by inflating the balloon, - Parietal fixation is carried out by two sutures thanks to the pre-positioned threads		3- Prosthesis in place: -The prosthesis is fixed, -The balloon must imperatively be withdrawn	

Description of the prosthesis:

This prosthesis is in a shape of a flat plug with 2 fixation threads crimped with needles (4 needles) for the product references PSUMBIN05D and PSUMBIN07D, and with 4 fixation threads crimped with needles (8 needles) for the PSUMBIN09D product reference.

The aim of its innovating expansion balloon concept is to make dissection easier and to secure the unfolding of the prosthesis.

This balloon is located inside the plug; it is filled-up with 30cm³ of air. It acts as a pneumatic Jack ensuring complete crease-free unfolding of the prosthesis.

Its double-sided design combines a smooth ePTFE side which reduces viscera adhesions with a knitted monofilament polypropylene side which allows better colonisation of the peritoneum.

Sizes of the prosthesis:

The medical device is a sterile kit.

The product range is composed of:

PSUMBIN 5cm diameter: for small hernias and ventral hernias on trocar site

PSUMBIN 7cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias < 2cm

PSUMBIN 9cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias that measure between 2 and 5cm

Installation of the patient :

Patient lying on his back

Both legs kept in the body's axis

Each arm placed on an arm wrester in a 90° angle from the body's axis

1-Incision :

The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia.

2-Dissection and resection of the peritoneal sac :

The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place.

Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand.

3-Positioning and fixating the prosthesis :

PROMESH® SURG UMB is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way.

The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice.

The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall.

Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall.

Be careful not to suture the balloon to the implant when fixing the prosthesis.

The balloon is deflated and removed.

4-Closure of the incision :

Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

R_x
only Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

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