**NOT FOR DISTRIBUTION IN THE
USA OR ITS TERRITORIES**

NOT291 Rev. C

Diam™ Intervertebral Prosthesis

2020-11-19

Medtronic

English

Diam™ Intervertebral Prosthesis

ARTICULAR PROCESS STABILIZATION AND DISCAL ASSISTANCE DEVICE

Note: not all parts may be available in each geography.

MATERIAL

- Dimethyl siloxane (silicone)
- Polyethylene terephthalate
- Titanium

INDICATIONS

Lumbar Surgery

- Arthropathic facet-syndromes.
- Posterior prosthesis for discal assistance.
- Transitional osteosynthesis.
- Foraminal stenosis.

The Diam™ Intervertebral Prosthesis allows an effective stabilization of the operated segments while ensuring a certain flexibility. The purpose of this stabilization is to decrease the lumbagos of the patient while avoiding the arthrodesis.

The Diam™ Intervertebral Prosthesis is designed for long-term implantation.

CONTRAINDICATIONS

Do not implant in case of intolerance to silicone, polyethylene terephthalate, or titanium. Do not implant during child growth.

Do not implant on septic zone or in case of septicemia.

POSSIBLE UNDESIRABLE SECONDARY EFFECTS

Known undesirable effects concerning interspinous devices unlikely related to the device, but related to the surgery or patient condition include:

- Inflammation.
- Permanent ligament injury
- Rupture of ligament.
- Removal of the prosthesis.
- Infection.
- Pseudo meningocele, fistula, persistent CSF leakage, and/or meningitis.
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anaesthesia, paraesthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits.

- Cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss.
- Urinary retention, loss of bladder control, or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery.
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: adjacent segment degeneration.
- Non-union or pseudarthrosis, delayed union.
- Mal-union.
- Cessation of any potential growth of the operated portion of the spine.
- Loss of or increase in spinal mobility or function.
- Inability to perform the activities of daily living. ▪ Bone loss or decrease in bone density.
- Graft donor site complications including pain, fracture, or wound healing problems.
- Ileus, gastritis, bowel obstruction, or loss of bowel control or other types of gastrointestinal system compromise.
- Haemorrhage, hematoma, occlusion, seroma, oedema, hypertension, embolism, stroke, excessive bleeding. ▪ Reproduction system compromise, sterility, sexual dysfunction.
- Development of respiratory problems (e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.). ▪ Change in mental status.
- Death.

Undesirable secondary effects related to the device include:

- Device loosening.
- Device breakage.
- Knitting tearing.
- Seroma/swelling around the device.
- Allergy.
- Foreign body reaction. ▪ Inflammation.
- Worsening pain.
- Inaccurate positioning (posterior or wrong level).

POSTOPERATIVE

When explanting and/or disposing of a device, avoid exposure to bodily substances such as blood, tissue, etc., as contact could lead to infection or disease. Always wear and use proper equipment, taking special care with sharp objects and needles. Follow your healthcare center's policy regarding both the disposal of devices and any events of exposure.

STORAGE

Keep in a dry place, without sunlight, and in an ambient atmosphere.

VISUAL INSPECTION

Visually inspect all sterile-barrier packaging before use. If the sterile barrier is damaged or the integrity has been compromised, do not use the product. Contact Medtronic for return information.

Visually inspect the device before use. If the device is damaged, do not use the product. Contact Medtronic for return information.

PACKAGING

Devices are supplied sterile. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package is broken, the product should not be re-sterilized. If a loaner set is used, all sets and components should be carefully checked for completeness and to ensure there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Medtronic.

WARNINGS

COUSIN SURGERY does not guarantee or recommend any special trademark for fixation devices. Properties of these devices are subjected to modifications made by the manufacturer and on which COUSIN SURGERY has no control.

Diam™ Intervertebral Prosthesis is not recommended for implantation in children. This product must be implanted only by a trained physician.

Physician note: although the physician is the learned intermediary between the company and the patient, the important medical information in this document should be conveyed to the patient.

This device was designed for single patient use only. Do not reprocess or reuse this product. Reuse or reprocessing may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

MRI INFORMATION

MR Conditional

Diam™ Intervertebral Prosthesis is MR Conditional.

The Diam™ Intervertebral Prosthesis was determined to be MR Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM

International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Product comparison demonstrated that the Diam™ Intervertebral Prosthesis is MR Conditional. A patient with this device can be scanned immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less.
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less.

MRI-Related Heating

MRI testing on Diam™ Intervertebral Prosthesis yields a similar result as compared to testing done on a worst case X-Stop implant during a 15 minute MRI scan using the 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI). The highest temperature change for the MR system was +1.7°C.

Therefore, the estimate based on comparison to a worst case condition (X-STOP™ Implant) MRI-related heating for the Diam™ Intervertebral Prosthesis at 3 Tesla using a transmit/receive RF body coil at an MR system reported whole body SAR averages of 2.9 W/kg (i.e. associated with a calorimetry measured whole body averaged value of 2.7 W/kg) and the greatest amount of expected heating that would occur in association with these specific conditions is equal to or less than +1.7°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the Diam™ Intervertebral Prosthesis. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

PRODUCT COMPLAINTS

To report product problems, contact Medtronic.

Patients in the European Union experiencing a serious incident in relation to the device should contact Medtronic and the competent authority of the Member State in which they are established.

FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is required, contact Medtronic.

©2020 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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Diam™ intervertebral protese er MR-betinget.

Diam™ intervertebral protese er funnet å være MR-betinget i henhold til terminologien som er angitt i American Society for Testing and Materials (ASTM) International, Designation: F2503-05, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

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