# PROMESH<sup>®</sup> SURG INTRA

## STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notices d'instructions</a>	Page	3
de	<a href="#">Gebrauchsanweisung</a>	Seite	4
it	<a href="#">Istruzioni per l'uso</a>	Pagina	5
es	<a href="#">Instrucciones de uso</a>	Página	6
pt	<a href="#">Nota de instruções</a>	Página	7
ελ	<a href="#">Οδηγίες χρήσης</a>	σελίδα	8
bg	<a href="#">Инструкции за употреба</a>	страни	9
da	<a href="#">Brugsvejledning</a>	Side	10
fi	<a href="#">Käyttöohjeet</a>	Sivu	11
hu	<a href="#">Felhasználási útmutató</a>	Pagina	12
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	13
pl	<a href="#">Instrukcja użycia</a>	Strona	14
ro	<a href="#">Návod k použití</a>	Pagină	15
sk	<a href="#">Návod k použiti</a>	Strana	16
cs	<a href="#">Pokyny k použití</a>	Strana	17
tr	<a href="#">Kullanım Kılavuzu</a>	Sayfa	18
ru	<a href="#">Инструкция по применению</a>	<u>страница</u>	19
zn	<a href="#">使用说明</a>	<u>页</u>	20
ar	<a href="#">تعليمات الاستخدام</a>	<u>صفحة</u>	21
sl	<a href="#">Navodila za uporabo</a>	<u>Stran</u>	22
sv	<a href="#">Bruksanvisning</a>	<u>Sida</u>	23
no	<a href="#">Bruksanvisning</a>	<u>Side</u>	24
uk	<a href="#">Інструкція із застосування</a>	<u>сторінка</u>	25

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**NOT282\_240920**

Version du 20/09/2024

This release is the last update of the instructions of use and replace the previous edition

## Instructions for use

en

### **IMPLANT PROMESH® SURG INTRA STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS SINGLE USE STERILE PRODUCT**

#### **DESCRIPTION**

PROMESH® SURG INTRA implants are parietal reinforcement implants.

#### **IMPLANTED MATERIALS**

Polypropylene and ePTFE for anti-adherent side (See on the package).

Biocompatible ink with dimethyl siloxane (depending on references) (See on the package).

**Origin neither human nor animal – Non resorbable.**

#### **INDICATIONS**

Repair and parietal reinforcement for umbilical and ventral hernias.

#### **PERFORMANCE**

These products can be used in laparoscopic surgery or open surgery. **PROMESH® SURG INTRA** products have a no-adhesive side, which is designed to avoid or reduce visceral adhesences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. Marking in black ink, composed of a cross and dots is on the non-adhesive side. The meshes, biocompatible and non resorbable have the advantage of having a shape memory and a great flexibility.

#### **CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

#### **UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- |                   |                      |                     |                            |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma          | - Foreign Body Reaction    |
| - Infection       | - Fistula formation  | - Mesh migration    | - Irritation near by organ |
| - Recurrence      | - Erosion            | - Mesh deformation  | - Obstruction              |
| - Inflammation    | - Seroma             | - Allergic reaction |                            |

#### **USAGE PRECAUTIONS**

**PROMESH® SURG INTRA** implants are delivered sterile. Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date.

COUSIN SURGERY does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

As indicated on the product labelling, the implant, is designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

#### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

Do not cut out the mesh and suture it at least 1cm of the edge. The trocar to be used must be chosen according to the size of the implant to be used. This is to allow passage without excessive force. Always roll the implant with the ePTFE face inside for the passage in the trocar. For use in keyhole surgery, we recommend the use of non-notched atraumatic forceps.

#### **STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

#### **IMPORTANT**

For more information about the use of this product, please contact your representative or distributor

Rx  
only

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician

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