# DIAM™

## Sterile Intervertebral Prostheses

---

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Made in France

**Rx**  
only

This release is the last update of the instructions of use and replace the previous edition  
**Product under regulation 2017/745 (EU).**

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### **IDENTIFICATION AND DESCRIPTION OF THE DEVICE**

The DIAM™ Prosthesis ensures a perfect intervertebral support between two adjoining spinous processes of the spine using a posterior surgical approach.

DIAM™ Prosthesis composed of a silicone wedge (dimethyl siloxane) which is surrounded by a knitted polyethylene terephthalate coating (PET) with two flat ligaments crimped with a needle at one end and a hardened spliced loop sewn at the other end. Each product contains also two crimps of diameter 6.5 mm made of titanium, the crimps are adapted to the flat ligaments.

The intervertebral Prosthesis available in 5 different sizes (8mm, 10mm, 12mm, 14mm and 16mm) of silicone wedge with the flat ligaments.

### **INDICATIONS**

DIAM™ prostheses are indicated in lumbar surgery in skeletally mature patients for:

- Arthropathic facet-syndrome
- Posterior prosthesis for discal assistance
- Transitional osteosynthesis
- Foraminal stenosis

The DIAM™ Intervertebral Prosthesis allows an effective stabilization of the operated segments while ensuring a certain flexibility. The purpose of this stabilization is to decrease the lumbagos of the patient while avoiding the arthrodesis.

### **CONTRAINDICATIONS**

- Do not implant in case of intolerance to silicone, polyethylene terephthalate, or titanium
- Do not implant during child growth
- Do not implant on septic zone or in case of septicaemia
- Do not implant in pregnant women
- Do not implant in patients with mental deficiencies that may alter their understanding of pre- and post-operative instructions.
- Do not implant in the presence of Severely distorted anatomy caused by congenital abnormalities or abdominal curvature of the backbone that do not allow the device to be implanted.
- Do not implant in the presence of a tumour process located in the implantation zone.
- Do not implant in case of recent trauma that interfere with the placement of the implant in the insertion zone.
- Do not implant in case of spondylolisthesis (all grades).
- Do not implant outside levels L3-L5.

Some factors may compromise successful implantation:

- Severe obesity (BMI greater than or equal to 40)
- Severe osteoporosis

### **UNDESIRABLE SIDE EFFECTS**

**1. Known undesirable effects concerning interspinous devices unlikely related to the device, but related to the surgery or patient condition include:**

- Inflammation

- Permanent ligament injury
- Rupture of ligament
- Removal of the prosthesis
- Infection
- Pseudo meningocele, fistula, persistent CSF leakage, and/or meningitis
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anaesthesia, paraesthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits
- Cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss
- Urinary retention, loss of bladder control, or other types of urological system compromise
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union
- Mal-union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Inability to perform the activities of daily living
- Bone loss or decrease in bone density
- Graft donor site complications including pain, fracture, or wound healing problems
- Ileus, gastritis, bowel obstruction, or loss of bowel control or other types of gastrointestinal system compromise
- Haemorrhage, hematoma, occlusion, seroma, oedema, hypertension, embolism, stroke, excessive bleeding.
- Reproduction system compromise, sterility, sexual dysfunction
- Development of respiratory problems (e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.).
- Change in mental status
- Death

**2. Undesirable secondary effects related to the device include:**

- Device loosening
- Device breakage
- Knitting tearing
- Seroma/swelling around the device
- Allergy
- Foreign body reaction
- Inflammation
- Worsening pain
- Inaccurate positioning (posterior or wrong level)

**TARGET POPULATION**

DIAM™ prostheses are recommended for implantation in adult (skeletally mature patients).

**IMPLANTED MATERIALS**

DIAM™ prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

<i>Reference Material</i>	CB001A0008 Size 8mm	CB001A0010 Size 10mm	CB001A0012 Size 12mm	CB001A0014 Size 14mm	CB001A0016 Size 16mm	CB001A0015*
Dimethyl siloxane (silicone)	3,53g ± 1,5%	3,84g ± 1,5%	4,16g ± 1,5%	4,47g ± 1,5%	4,78g ± 1,5%	0g
Polyethylene terephthalate	2,62g ± 1,5%	2,67g ± 1,5%	2,73g ± 1,5%	2,78g ± 1,5%	2,83g ± 1,5%	0g

Titanium	0,38g ± 5%	0,38g ± 5%	0,38g ± 5%	0,38g ± 5%	0,38g ± 5%	0,38g ± 5%
----------	------------	------------	------------	------------	------------	------------

\*6,5mm extra crimps for DIAM™ prosthesis with flats ligaments.

**Origin neither human nor animal – Non resorbable.**

**PERFORMANCES**

DIAM™ ensures a certain flexibility and effective stabilization of operate segments in both flexion and extensions. This stabilization is intended to reduce the patient's lumbagos while avoiding arthrodesis. The motion scores are claimed as good as preoperatively.

The use of DIAM™ in these indications is designed to relieve pressure and reduce lower back pain related to lumbar pathologies.

**LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS :**

[https://www.ifu-oem.com/d-6318526d8b749\\_under-construction.pdf](https://www.ifu-oem.com/d-6318526d8b749_under-construction.pdf)

**MRI SAFETY**

**MRI INFORMATION**

**Conditional**

In non-clinical testing the DIAM™ **Spinal Stabilization System** was determined to be MRConditional. A patient with this device can be safely scanned immediately after device placement under the following conditions:

**Static Magnetic Field**

- Static magnetic field of 1.5 Tesla and 3-Tesla.
- Maximum spatial gradient magnetic field of 3000-Gauss/cm or less
- Maximum whole body average specific absorption rate (SAR) of 4.0 -W/kg or less under normal mode or first level controlled operating mode, for 15 minutes of scanning per pulse sequence.

**MRI-Related Heating**

Based on measurements of RF heating according to ASTM F2182-11a, the DIAM™ Spinal Stabilization System is expected to produce a maximum temperature of 3.2°C for a whole body SAR of 4.0 W/kg in a 1.5-Tesla/64-MHz and 3.0-Tesla/128-MHz MRI Systems, for a 15-min of scanning (i.e., per pulse sequence)."

**Artifact Information**

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the DIAM™ Spinal Stabilization System. The maximum artifact beyond the implant was 5 mm for the spin echo sequence and 9 mm for the gradient echo sequence in a 3-Tesla MR system (GE Signa HDxt MR System). Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

The presence of other implants or the health state of the patient may require a modification of the MR conditions.

**MRI Patient Counseling Information**

Physicians should communicate with the patient the following information about Magnetic Resonance Imaging (MRI) with respect to DIAM™ **Spinal Stabilization System**:

DIAM™ **Spinal Stabilization System** performance has been established for MRI systems at field strengths of 1.5 Tesla and 3.0 Tesla only. During and MRI the patient may notice a warming sensation around the implant or feel a tingling sensation. If the warming or tingling sensation is uncomfortable the patient should communicate this to the MR technologist, the MRI should be stopped and the settings adjusted to reduce or eliminate the sensation. The highest temperature change above the background observed in non-clinical testing was +3.2°C (associated with specific conditions listed above). Additionally, the metal in the implant may cause the MRI image to be distorted in the area around the implant. The MRI can be adjusted to minimize the image distortion.

Physicians should instruct patients to:

- Inform any healthcare personnel (e.g. doctor or MR technologist) that they have an implanted DIAM™ Spinal Stabilization System prior to receiving an MRI.
- The patient's doctor will recommend whether or not an MRI is appropriate.

**PRECAUTIONS FOR USE**

The DIAM™ devices are delivered sterile (gamma rays sterilization).

Before any use, inspect the integrity of the packaging and device (including pouches).

Do not use in the event of deterioration of the labels and/or the device and/or the packaging.  
Do not use if the device presents malfunction, please use another device, and refers to the section INFORMATION REQUEST AND COMPLAINTS of this notice in order to notify it.  
Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.  
When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.  
Do not use if the device is out of date.  
DIAM™ implant is intended to be only implanted with specific instruments. These are supplied and marked CE by Companion Spine.

The DIAM™ prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of DIAM™. The information contained in this notice is necessary but not sufficient to control the surgical technique.

The prescription of the device is decided by the surgeon, the only person qualified to do so.

The DIAM™ is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

Topping off is not contra-indicated but special attention is required. It is not recommended to place more than one DIAM™ implant; The DIAM™ must not be placed in contact with another device on the same level. It is possible to combine the DIAM™ with another device but not on the same level to preserve integrity. Any contact could damage the implant. Once the needle is through the interspinous ligament, grab directly the tethers using a suitable suturing instrument (needle holder/tweezers) to slide I completely through the ligament to then apply tension. Do not use the needle to pull or apply tension.

The biocompatibility assessment has been performed and validated for a permanent implantation.

The mechanical, clinical and PMS elements guarantee the performance of the product for a duration of 5 years.

#### **IMPORTANT: DO NOT REUSE AND DO NOT RESTERILIZE THE DIAM™ PROSTHESIS**

As indicated on the product labelling, the prosthesis is designed for single use. It cannot be re-used and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

#### **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The patient should be informed of the post-operative care and precautions to be observed.

The use of an orthopaedic appliance may be necessary but is not mandatory. This decision can only be made by the surgeon.

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal.

An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. Regarding the lifetime of the product once implanted indicated on the implant card, Cousin Biotech recommends that the surgeon explain to the patient that the implant is designed for permanent implantation. The duration indicated on the card is indicative based on mechanical and clinical data. Please refer to "the Precaution for use" of this instruction for use about the lifetime. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

#### **STORAGE PRECAUTIONS**

Keep in a dry place away from sunlight, at room temperature, in its original packaging.

Cousin Biotech carried out tests on the products, which were placed at temperatures and humidity levels ranging from 18 to 40 degrees and from 50% to 90% of humidity. The results were compliant, with no impact on the products.

#### **LINK TO SURGICAL TECHNIQUE**

The surgical technical is available in electronic format on the website: [https://www.ifu-oem.com/d-6318526d8b749\\_under-construction.pdf](https://www.ifu-oem.com/d-6318526d8b749_under-construction.pdf)

#### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any

implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

Risks associated with the explanation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explanation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

#### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

COUSIN BIOTECH does not guarantee or recommend any special trademark for fixation devices. The properties of these devices are subject to change by the manufacturer over which COUSIN BIOTECH has no control.

#### **INFORMATION REQUEST AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to produce every effort and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question and inform COUSIN BIOTECH without delay. For any correspondence, please specify the reference, batch number, reference contact details and a comprehensive description of the incident or complaint.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor. Brochures, documentation, and surgical technique are available on request from COUSIN BIOTECH and Companion Spine. To report product problems, contact Companion Spine.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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